


Product Datasheet

Human Anti-HK1(anti-Hexokinase-1 antibody)ELISA Kit (orb1715126)

Description	Human Anti-HK1(anti-Hexokinase-1 antibody)ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	Anti-HK1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P19367
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
Anti-
HK1(anti-
Hexokinase-
1 antibod...