


Product Datasheet

Human Nectin 2(Nectin cell adhesion molecule 2) ELISA Kit (orb1715131)

Description	Human Nectin 2(Nectin cell adhesion molecule 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	Nectin 2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q92692
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
Nectin
2(Nectin
cell
adhesion
mole...