


## Product Datasheet

### Mouse GHRL(Ghrelin) ELISA Kit (orb1715132)

|                        |                               |
|------------------------|-------------------------------|
| <b>Description</b>     | Mouse GHRL(Ghrelin) ELISA Kit |
| <b>Reactivity</b>      | Mouse                         |
| <b>Range</b>           | 3.125-200pg/ml                |
| <b>Target</b>          | Ghrelin                       |
| <b>Storage</b>         | 2-8 °C for 6 months           |
| <b>Note</b>            | For research use only         |
| <b>Uniprot ID</b>      | <b>Q9EQX0</b>                 |
| <b>Sensitivity</b>     | 1.875pg/ml                    |
| <b>Expiration Date</b> | Please enquire.               |



Mouse  
GHRL(Ghrelin)  
ELISA Kit