


Product Datasheet

Mouse Nqo1(NAD(P)H dehydrogenase [quinone] 1) ELISA Kit (orb1715133)

Description	Mouse Nqo1(NAD(P)H dehydrogenase [quinone] 1) ELISA Kit
Reactivity	Mouse
Range	46.875-3000pg/ml
Target	Nqo1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q64669
Sensitivity	28.125pg/ml
Expiration Date	Please enquire.



Mouse
Nqo1(NAD(P)H
dehydrogenase
[quinon...