


Product Datasheet

Rat Sox9(Transcription factor SOX-9) ELISA Kit (orb1715141)

Description	Rat Sox9(Transcription factor SOX-9) ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Sox9
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	F1LYL9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Sox9(Transcription
factor SOX-9)
ELI...