


Product Datasheet

Human QSOX1(Sulfhydryl oxidase 1) ELISA Kit (orb1715148)

Description	Human QSOX1(Sulfhydryl oxidase 1) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	QSOX1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O00391
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
QSOX1(Sulfhydryl
oxidase 1) ELISA
...