


Product Datasheet

Human GPM6A(Neural membrane glycoprotein M6-a) ELISA Kit (orb1715153)

Description	Human GPM6A(Neural membrane glycoprotein M6-a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	GPM6A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P51674
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
GPM6A(Neural
membrane
glycoprote...