


Product Datasheet

Human CHRNA1(Acetylcholine receptor subunit alpha)ELISA Kit (orb1715157)

Description	Human CHRNA1(Acetylcholine receptor subunit alpha)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	CHRNA1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P02708
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
CHRNA1(Acetylcholine
receptor subu...