


Product Datasheet

Human YTHDF3(YTH domain-containing family protein 3) ELISA Kit (orb1715167)

Description	Human YTHDF3(YTH domain-containing family protein 3) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	YTHDF3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7Z739
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
YTHDF3(YTH
domain-
containing
famil...