


Product Datasheet

Mouse Trpv1(Transient receptor potential cation channel subfamily V member 1) ELISA Kit (orb1715168)

Description	Mouse Trpv1(Transient receptor potential cation channel subfamily V member 1) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Trpv1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q704Y3
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Trpv1(Transient
receptor
potential...