


Product Datasheet

Human ATR(Serine/threonine-protein kinase ATR) ELISA Kit (orb1715170)

Description	Human ATR(Serine/threonine-protein kinase ATR) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	ATR
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13535
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
ATR(Serine/threonine-
protein kinas...