


Product Datasheet

Human EIF2C1(Protein argonaute-1) ELISA Kit (orb1715177)

Description	Human EIF2C1(Protein argonaute-1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	EIF2C1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UL18
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
EIF2C1(Protein
argonaute-1)
ELISA ...