


Product Datasheet

Human CCL26(C-C motif chemokine 26) ELISA Kit (orb1715183)

Description	Human CCL26(C-C motif chemokine 26) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	CCL26
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y258
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
CCL26(C-C
motif
chemokine
26) ELIS...