


Product Datasheet

Human URGCP(Up-regulator of cell proliferation) ELISA Kit (orb1715186)

Description	Human URGCP(Up-regulator of cell proliferation) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	URGCP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8TCY9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
URGCP(Up-
regulator
of cell
prolife...