


Product Datasheet

Mouse Camk2b(Calcium/calmodulin-dependent protein kinase type II subunit beta) ELISA Kit (orb1715189)

Description	Mouse Camk2b(Calcium/calmodulin-dependent protein kinase type II subunit beta) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Camk2b
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P28652
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Camk2b(Calcium/calmodulin-
dependen...