


Product Datasheet

Rabbit PPIA(Cyclophilin A) ELISA Kit (orb1715192)

Description	Rabbit PPIA(Cyclophilin A) ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/ml
Target	PPIA
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9TTC6
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rabbit
PPIA(Cyclophilin
A) ELISA Kit