


Product Datasheet

Mouse Gpt(Alanine aminotransferase 1) ELISA Kit (orb1715193)

Description	Mouse Gpt(Alanine aminotransferase 1) ELISA Kit
Reactivity	Mouse
Range	0.313-20mIU/ml
Target	Gpt
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8QZR5
Sensitivity	0.188mIU/ml
Expiration Date	Please enquire.



Mouse
Gpt(Alanine
aminotransferase
1) EL...