


Product Datasheet

Porcine ADP/Acrp30(Adiponectin) ELISA Kit (orb1715196)

Description	Porcine ADP/Acrp30(Adiponectin) ELISA Kit
Reactivity	Porcine
Range	3.125-200pg/ml
Target	ADP/Acrp30
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q6PP07
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.



Porcine
ADP/Acrp30(Adiponectin)
ELISA Ki...