


## Product Datasheet

### Human BRN2(Brain-specific homeobox/POU domain protein 2)ELISA Kit (orb1715198)

|                        |                                                                   |
|------------------------|-------------------------------------------------------------------|
| <b>Description</b>     | Human BRN2(Brain-specific homeobox/POU domain protein 2)ELISA Kit |
| <b>Reactivity</b>      | Human                                                             |
| <b>Range</b>           | 15.625-1000pg/ml                                                  |
| <b>Target</b>          | BRN2                                                              |
| <b>Storage</b>         | 2-8 °C for 6 months                                               |
| <b>Note</b>            | For research use only                                             |
| <b>Uniprot ID</b>      | <b>P20265</b>                                                     |
| <b>Sensitivity</b>     | 9.375pg/ml                                                        |
| <b>Expiration Date</b> | Please enquire.                                                   |



Human  
BRN2(Brain-  
specific  
homeobox/POU  
d...