


Product Datasheet

Bovine Beta-LG(beta-lactoglobulin)ELISA kit (orb1715200)

Description	Bovine Beta-LG(beta-lactoglobulin)ELISA kit
Reactivity	Bovine
Range	0.156ug/ml-10ug/ml
Target	Beta-LG
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P02754
Sensitivity	0.094ug/ml
Expiration Date	Please enquire.



Bovine Beta-LG(beta-lactoglobulin) ELISA...