

Product Datasheet

anti- TMEM173 antibody (orb1715205)

Description	anti- TMEM173 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	transmembrane protein 173
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	39 kDa
Uniprot ID	Q86WV6
Dilution Range	WB: 1:1000 - 1:5000; IHC: 1:100 - 1:500
Expiration Date	12 months from date of receipt.