

Product Datasheet

B7-H3 Antibody / CD276 (orb1716146)

Description	T cell activation and immune function are regulated by the innate immune system through positive and...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 100-300 from the human CD276 protein was used as the immunogen for the B7-H3 antibody.
Storage	Aliquot the B7-H3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the B7-H3 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q5ZPR3
Hazard Information	This B7-H3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.