

Product Datasheet

PDCD1 Antibody / PD-1 / PD1 (orb1716185)

Description	PDCD-1 (programmed cell death-1 protein), also designated CD279, is a type I transmembrane receptor ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A portion of amino acids 1-200 from the human protein was used as the immunogen for the recombinant PDCD1 antibody.
Storage	Aliquot the recombinant PDCD1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant PDCD1 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q15116
Hazard Information	This recombinant PDCD1 antibody is available for research use only.
Dilution Range	ELISA: for coating, order antibody without BSA, Immunohistochemistry (FFPE): 1-2 µg/ml
Expiration Date	12 months from date of receipt.