

Product Datasheet

Golgi Antibody / GLG1 (orb1716243)

Description

The Golgi apparatus is an organelle present

endomem...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC

Applications

Immunogen

Recombinant full-length human protein was used as the immunogen for the Golgi Glycoprotein 1 antibody.

Storage

Aliquot the Golgi Glycoprotein 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Optimal dilution of the Golgi Glycoprotein 1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID

Q92896

Hazard

Information

This Golgi Glycoprotein 1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (Frozen & FFPE): 0.5-1ug/ml (No special pretreatment is required for the immunohistochemical staining of formalin-fixed tissues)

Expiration Date

12 months from date of receipt.