

Product Datasheet

ATRX Antibody / RAD54 (orb1716244)

Description	ATRX is a member of the Snf2 family of helicase/ATPases, which contribute to the remodeling of the n...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 2200-2450 from the human protein was used as the immunogen for the ATRX antibody.
Storage	Aliquot the ATRX antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the ATRX antibody should be determined by the researcher.
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P46100
Hazard Information	This ATRX antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.