

Product Datasheet

Steroidogenic Factor 1 Antibody / SF-1 / NR5A1 (orb1716257)

Description

Steroidogenic Factor 1 (SF-1) is

that belongs to subfamily 5. ...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A portion of amino acids 220-461 from the human protein was used as the immunogen for the recombinant Steroidogenic Factor 1 antibody.

Storage

Aliquot the recombinant Steroidogenic Factor 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Optimal dilution of the recombinant Steroidogenic Factor 1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Uniprot ID

[Q13285](#)

Hazard
Information

This recombinant Steroidogenic Factor 1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.