

Product Datasheet

TRPV1 polyclonal antibody (orb1716327)

Description	TRPV1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human TRPV1(NP_061197.4).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q8NER1
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.