

Product Datasheet

Cleaved PARP (Asp214) polyclonal antibody (orb1716475)

Description	Cleaved PARP (Asp214) polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of Mouse Cleaved PARP(Asp214).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	Q921K2
Dilution Range	WB,1:1000 - 1:2000
Expiration Date	12 months from date of receipt.