

Product Datasheet

Cleaved PARP (Asp214) polyclonal antibody (orb1716476)

Description	Cleaved PARP (Asp214) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Cleaved PARP (Asp214)(Asp214).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	89KDa
Uniprot ID	P09874
Dilution Range	WB, 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.