

Product Datasheet

BrdU polyclonal antibody (orb1716606)

Description	BrdU polyclonal antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	BrdU
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, 1:50 - 1:500 IF/ICC, 1:50 - 1:500
Expiration Date	12 months from date of receipt.