

Product Datasheet

SARS-CoV-2 Spike S1 polyclonal antibody (orb1716618)

Description	SARS-CoV-2 Spike S1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	DOT, IF, IP, WB
Immunogen	Recombinant fusion protein coronavirus Spike S1(YP_009724390.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	110KDa
Uniprot ID	P0DTC2
Dilution Range	DB,1:500 - 1:2000 ELISA,1:20000-1:80000 WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 IP,1:50 - 1:200
Expiration Date	12 months from date of receipt.