

Product Datasheet

SMAD2/SMAD3 polyclonal antibody (orb1716640)

Description	SMAD2/SMAD3 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated of human SMAD2/SMAD3SMAD3(NP_005893.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	52kDa/60kDa
Uniprot ID	P84022, Q15796
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.