

Product Datasheet

LAMP1 polyclonal antibody (orb1716766)

Description	LAMP1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human LAMP1(NP_005552.3).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	45KDa/110KDa
Uniprot ID	P11279
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.