


Product Datasheet

RNF17 antibody (orb1717)

Description	RNF17 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RNF17 (1031-1200/1623 aa)
Target	RNF17
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	185 kDa
Uniprot ID	Q9BXT8
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat test...