

Product Datasheet

MLKL polyclonal antibody (orb1717196)

Description	MLKL polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Recombinant fusion protein of mouse MLKL(Q9D2Y4).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	Q8NB16
Dilution Range	IHC,1:50 - 1:200 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.