

Product Datasheet

LC3B polyclonal antibody (orb1717223)

Description	LC3B polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human LC3B(NP_073729.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	14KDa/16KDa
Uniprot ID	Q9GZQ8
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.