

Product Datasheet

VDAC1 / Porin polyclonal antibody (orb1717287)

Description	VDAC1 / Porin polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human VDAC1 / Porin(NP_003365.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	32KDa
Uniprot ID	P21796
Dilution Range	WB,1:500 - 1:2000 IHC,1:100 - 1:200
Expiration Date	12 months from date of receipt.