

Product Datasheet

Phospho-p90Rsk/RSK1/RPS6KA1-T573 polyclonal antibody (orb1717379)

Description	Phospho-p90Rsk/RSK1/RPS6KA1-T573 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T573 of human p90Rsk/RSK1/RPS6KA1(NP_002944.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	Q15418
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.