

Product Datasheet

Phospho-Na⁺/K⁺-ATPase-S16 polyclonal antibody (orb1717389)

Description	Phospho-Na ⁺ /K ⁺ -ATPase-S16 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S16 of human Na ⁺ /K ⁺ -ATPase(NP_000692.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	112kDa
Uniprot ID	P05023
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.