

Product Datasheet

Phospho-Histone H3-S10 polyclonal antibody (orb1717416)

Description	Phospho-Histone H3-S10 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around S10 of human Histone H3(NP_003484.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	17KDa
Uniprot ID	Q16695
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.