

Product Datasheet

Phospho-POLR2A CTD-S5 polyclonal antibody (orb1717420)

Description	Phospho-POLR2A CTD-S5 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S5 of human POLR2A CTD repeatYSPTSPS.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	270kDa
Uniprot ID	P24928
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:100
Expiration Date	12 months from date of receipt.