

Product Datasheet

Phospho-p70 S6 Kinase 1-T389 polyclonal antibody (orb1717457)

Description	Phospho-p70 S6 Kinase 1-T389 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic phosphorylated peptide around T389 of human P70 S6K(NP_003152.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	70KDa
Uniprot ID	P23443
Dilution Range	WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:200
Expiration Date	12 months from date of receipt.