

Product Datasheet

Phospho-p90Rsk/RSK1/RPS6KA1-S380 polyclonal antibody (orb1717458)

Description	Phospho-p90Rsk/RSK1/RPS6KA1-S380 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S380 of human p90Rsk/RSK1/RPS6KA1(NP_002944.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	90kDa
Uniprot ID	Q15418
Dilution Range	WB,1:500 - 1:2000 IP,1:50 - 1:100
Expiration Date	12 months from date of receipt.