

Product Datasheet

Phospho-STAT3-S727 polyclonal antibody (orb1717483)

Description	Phospho-STAT3-S727 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic phosphorylated peptide around S727 of human Phospho-STAT3-S727(NP_644805.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	86KDa
Uniprot ID	P40763
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200 IP,1:50 - 1:100
Expiration Date	12 months from date of receipt.