

Product Datasheet

Phospho-ERK1-T202 + ERK2-T185 polyclonal antibody (orb1717502)

Description	Phospho-ERK1-T202 + ERK2-T185 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around T202 of human ERK1/ERK2(NP_002737.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	42KDa/44KDa
Uniprot ID	P27361, P28482
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:100 IF/ICC,1:20 - 1:100
Expiration Date	12 months from date of receipt.