

Product Datasheet

Phospho-PLC gamma 1 (PLCG1)-Y783 polyclonal antibody (orb1717508)

Description	Phospho-PLC gamma 1 (PLCG1)-Y783 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y783 of human PLC gamma 1 (PLCG1)(NP_877963.1?).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	155KDa
Uniprot ID	P19174
Dilution Range	WB, 1:500 - 1:2000
Expiration Date	12 months from date of receipt.