

Product Datasheet

Phospho-STAT1-Y701 polyclonal antibody (orb1717510)

Description	Phospho-STAT1-Y701 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y701 of human STAT1(NP_009330.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	84KDa/91KDa
Uniprot ID	P42224
Dilution Range	WB,1:500 - 1:1000 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.