

Product Datasheet

Acetyl-Histone H4-K5/K8/K12/K16 polyclonal antibody (orb1717526)

Description	Acetyl-Histone H4-K5/K8/K12/K16 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic acetylated peptide around K5/K8/K12/K16 of human HistoneH4.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	P62805
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.