

Product Datasheet

ROBO3 polyclonal antibody (orb1718259)

Description	ROBO3 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ROBO3(NP_071765.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	200kDa
Uniprot ID	Q96MS0
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.