

Product Datasheet

FOXC1 polyclonal antibody (orb1718457)

Description	FOXC1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FOXC1(NP_001444.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q12948
Dilution Range	WB, 1:100 - 1:1000
Expiration Date	12 months from date of receipt.